

SBS Format Combo-NEW

I. APPLICATIONS

1. Sample collection
2. Sample transportation
3. Sample storage
4. Dispensing samples

II. FEATURES

1. Racks: 2D barcode and human readable numbers on the side, 2D barcode on the bottom; Vials: 2D barcode on the bottom
2. Simple and powerful mechanical locking performance, slide the buttons on two sides of the rack to unlock
3. Blue capped vials; white racks, 96 square wells
4. Working temperature: stable from -196°C to 121°C
5. Suitable for liquid nitrogen freezing

Warning:

Do not use cryogenic vials in liquid phase of liquid nitrogen. Ideal for cryogenic sample storage in gas phase of liquid nitrogen.

III. DESIGN

1. Vials, caps and a rack in one package, convenient to use
2. Fours types of codes on racks and 2D barcodes and human-readable codes on vials enable quick scan and easy identification
3. Vacuum packaging, foreign DNA free, DNase & RNase free, gamma radiation.

IV. PACKAGING

Item	Description	Packaging	Sterile
89-4233	0.3ml Vials+Rack, 96-well	1 Sets/Bag , 10 Bags/Pack , 2 Packs/Case	Yes
89-4253	0.5ml Vials+Rack, 96-well	1 Sets/Bag , 10 Bags/Pack , 2 Packs/Case	Yes
89-4283	0.75ml Vials+Rack, 96-well	1 Sets/Bag , 10 Bags/Pack , 2 Packs/Case	Yes
89-4213	1ml Vials+Rack, 96-well	1 Sets/Bag , 10 Bags/Pack , 2 Packs/Case	Yes

V. STERILITY

For each lot, sterilization records are recorded and signed off by qualified personnel for product release.

VI. CERTIFICATION




Biologix products are manufactured completely according to the TUV ISO 9001:2015 standard.

VII. PRODUCT PHOTOS




VIII. PRODUCT PARAMETER


Dimensions=mm



Caps of 89-4233 & 89-4253



89-4233 & 89-4253



89-4283 & 89-4213