

## Drosophila Bottle

### I. APPLICATIONS

1. Storage accessory for Drosophila related research

### II. FEATURES

1. Materials: polypropylene
2. Autoclavable
3. Disposable
4. Transparent for easy sample viewing

### III. BOTTLE DESIGN

1. 6oz (177ml) bottle
2. Molded graduations in 25ml increments up to 100ml
3. Square bottom
4. Available in both bulk and tray packaging

### IV. PACKAGING

1. Bulk, 200 pieces/case
2. Tray, 25 bottles/tray, 200 pieces/case

### V. DESCRIPTION

Item	Description	Packaging	Sterile
51-17720	Drosophila bottles, square bottom	Bulk, 200 pieces/case	No
51-17721	Drosophila bottles, square bottom	Tray, 25/tray, 8 trays/case	No

### VI. STERILITY

1. N/A

### VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.


### VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

### IX. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G	H
51-17720	Drosophila Bottle	103.3	74.6	12.0	6.0	O.D. 43.0	O.D.38.5	55.6	34.2



51-17720