

Cell Spreader

I. APPLICATIONS

1. Microbiology
2. Spreading of cells across growth surfaces

II. FEATURES

1. Materials: polystyrene
2. Gamma radiation sterilized
3. Disposable

III. SPREADER DESIGN

1. L-Shaped design
2. Available individually wrapped or 10 pieces/pack
3. Smooth surface for even spreading

IV. PACKAGING

1. Bulk, Individually wrapped, 500 pieces/case
2. Bulk, 10 pieces/pack, 50 packs/case

V. DESCRIPTION

Item	Description	Packaging	Sterile
65-1001	Cell spreader, orange	Individually wrapped, 500 pieces/case	Yes
65-1010	Cell spreader, orange	10 pieces/pack, 500 pieces/case	Yes

VI. STERILITY

1. Products are Gamma Radiation Sterilized according to ISO13485 standard.
2. For each lot, sterilization records are recorded and signed off by qualified personnel for product release.

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C	D	E	F
65-1001	Cell Spreader	144.0	33.2	O.D.4.0	O.D. 5.3	135°	O.D. 4.0

