

## CryoKING® PC Cryogenic Boxes

### I. APPLICATIONS

1. Storage and organization of microtubes/cryogenic vials

### II. FEATURES

1. Materials: durable polycarbonate
2. Designed for all CryoKing® cryogenic vials
3. Scanner-ready design for CryoKing® 2D barcode scanner and compatible with all other cryobox scanners
4. Scanner-ready design for CryoKing® RFID scanners
5. Chemically resistant to alcohols and mild organic solvents
6. Temperature range: stable from -196°C to 121°C (Autoclavable)
7. Designed & tested for long-term liquid nitrogen storage & Deep Freezer Use
8. Always 10 x 10 design to store 23% more samples compared with 9 x 9 boxes

### III. BOX DESIGN

1. 10x10 well storage (100-wells total)
2. The hollow bottom design, make the CryoKING® Cryogenic Vials be scanned easily
3. 1-inch and 2-inch Designs Available
4. Blue Color box with transparent lid
5. Holds 0.5ml, 1.0ml & 1.5ml microtubes/ cryogenic vials
6. Removable lid includes a numbered grid for orientation

### IV. PACKAGING

1. Bulk pack, 3 Pieces/Bag, 2 Bags/Pack, 2 Packs/Case, 12 Boxes/Case

### V. DESCRIPTION

Item	Description	Packaging	Sterile
98-0113	1-Inch Cryogenic storage box, Polycarbonate, 100-well, Blue	3 Pieces/Bag, 2 Bags/Pack, 2 Packs/Case	No
98-0213	2-Inch Cryogenic storage box, Polycarbonate, 100-well, Blue	3 Pieces/Bag, 2 Bags/Pack, 2 Packs/Case	No

### VI. STERILITY

1. N/A

**VII. CERTIFICATION**

Biologix products are manufactured completely according to the TUV ISO 9001:2015 standard.

**VIII. RESTRICTIONS**

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.


**IX. PRODUCT PHOTOS**



**98-0113**




**98-0213**




**X. PRODUCT PARAMETER**


Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G	H	I
98-0113	Box	21.75	15.5	5	129	130.1	131.8	11.3	132.5	19.7
98-0213	Box	21.75	15.5	5	129	130.1	131.8	11.3	132.5	35.9



98-0113





H

OCERVO® KING		3	4	5	6	7	8	9	10
1	2								
11	12	13	14	15	16	17	18	19	20
21	22	23	24	25	26	27	28	29	30
31	32	33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48	49	50
51	52	53	54	55	56	57	58	59	60
61	62	63	64	65	66	67	68	69	70
71	72	73	74	75	76	77	78	79	80
81	82	83	84	85	86	87	88	89	90
91	92	93	94	95	96	97	98	99	100



98-0213