

Product Name: Low-adhesion pipette tips & Tip Boxes	
Product Application:	Pipette tips with in boxes Googlab Professional Line are ultra-clear and fit universally the most pipettes. Ideal for transferring or measuring out liquid during experiment, especially in a laboratory.
Product Materials:	Polypropylene (PP). Meets USP, Class VI standards
Product Specifications:	Temperature range: stored at room temperature Shelf life: 3 Years after date of production (ensure package is in good)
Sterilization:	Yes (for Boxed tips). Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
RNase/DNase: Pyrogens:	Non-Pyrogenic DNase/RNase free
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Leak test: qualified Pipette Fit Test: qualified Packaging inspection: qualified
Features:	<ol style="list-style-type: none"> 1. Both tips and boxes are autoclavable. 2. With Graduation. 3. With precision mold and professional injection skill, Low-adhesion tips in boxes Googlab Professional Line has smooth inner surface to ensure Low liquid attachment and accurate sample results during liquid transfer. 4. Resistant to organic solvents. 5. Ultra-clear polypropylene good for observation. 6. Boxed tips are DNase/RNase free, and pyrogenic free. 7. Low-adhesion tips are compatible with most of the popular pipettes on the market. 8. Compliant with the USP 3127 standard, i.e. free from heavy metals. 9. Safe endotoxin level, compliant with USP/FDA standards High clarity and compatibility. 10. Perfect work with multi-channel pipettes. 11. High tightness on the pipette. 12. Resistance to organic solvents. 13. High protection against cross contamination.

Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland

E-mail: office@genoplast.com

Warehouse in the USA

GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com

Features:

- 14. Full production automation and Clean Room conditions.
- 15. Bagged tips are non-sterilized.
- 16. Boxed tips are sterilized by E-beam.


Product Range:

Low-adhesion pipette tips & Tip Boxes


SKU	Volume (µl)	Sterile	Material	LENGTH	Type of Packaging	Pack	Case
GBPT0010-B-N-LB	10µl	Yes	PP	40mm	Bagged	1x1000pcs.	1000pcs.
GBPT0010-R-NS	10µl	Yes	PP	40mm	Box (96pcs.)	10x96pcs.	960pcs.
GBPT0200-B-Y-LB	200µl	Yes	PP	52mm	Bagged	1x1000pcs.	1000pcs.
GBPT0200-R-NS	200µl	Yes	PP	52mm	Box (96pcs.)	10x96pcs.	960pcs.
GBPT1000-B-N-LB	1000µl	Yes	PP	88mm	Bagged	1x1000pcs.	1000pcs.
GBPT1000-R-NS	1000µl	Yes	PP	88mm	Box (96pcs.)	10x96pcs.	960pcs.

Technical Drawing:


10µl



200µl



1000µl



Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland

E-mail: office@genoplast.com

Warehouse in the USA

GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com