

Product Name:	Petri Dish
Product Application:	Suitable for bacterial culture.
Product Materials:	Polystyrene (PS)
Product Specifications:	Temperature range: stored at room temperature Shelf life: 5 Years after date of production (ensure package is in good)
Sterilization:	Yes. Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
Pyrogens:	Non-Pyrogenic
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: pass Packaging inspection: pass
Features:	<ol style="list-style-type: none"> 1. High clarity polystyrene 2. Flat transparent surface 3. Numeric mark on the bottom for easy identification 4. PA/PE packaging 5. Vent for gas exchange 6. Stackable for easy storage and handling 7. Non-pyrogenic 8. Sterilized by E-Beam, SAL=10⁻⁶

Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland

E-mail: office@genoplast.com

Warehouse in the USA

GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com


Product Range:

Petri Dish

SKU	Material	Diameter	Height	Ventilation ribs	Sterile	Pack	Case
GP521	PS	140mm	20mm	No	Yes	11x10pcs.	110pcs.
GP531	PS	140mm	20mm	Yes	Yes	11x10pcs.	110pcs.
GP506	PS	90mm	14,2mm	No	Yes	24x25pcs.	600pcs.
GP508	PS	90mm	14,2mm	Yes	Yes	24x25pcs.	600pcs.
GP507	PS	90mm	16,2mm	No	Yes	24x25pcs.	600pcs.
GP509	PS	90mm	16,2mm	Yes	Yes	24x25pcs.	600pcs.
GP516	PS	60mm	14,2mm	No	Yes	67x15pcs.	1005pcs.
GP504	PS	60mm	14,2mm	Yes	Yes	67x15pcs.	1005pcs.

Technical Drawing:

GP521 140mm



Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland



E-mail: office@genoplast.com

Warehouse in the USA

GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com

GP531 140mm



Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland


E-mail: office@genoplast.com

Warehouse in the USA



GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com

GP506 90mm



GP506 90mm



Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland


E-mail: office@genoplast.com

Warehouse in the USA



GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com

GP507 90mm



GP509 90mm



Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland


E-mail: office@genoplast.com

Warehouse in the USA



GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com

GP506 90mm



GP506 90mm



Warehouse in Europe

GenoPlast Biotech S.A.
Street. Brzozowa 8
83-200 Rokocin
Poland

E-mail: office@genoplast.com

Warehouse in the USA

GoogLab Scientific USA LLC
969 2nd Street SE
Charlottesville, VA 22902

E-mail: office@googlab.com